

Virtual environment

Installation

If you are using usegalaxy, conda is already installed. Skip the installation and continue with practice.

1. Get the conda installer

```
#Find the miniconda version you want to install in https://docs.conda.io/en/latest/miniconda.html, right click and copy link.  
#Use "wget" + link to download file  
#Example  
mkdir program  
cd program  
wget https://repo.anaconda.com/miniconda/Miniconda3-py38_23.5.2-0-Linux-aarch64.sh
```

2. Run the installer

```
bash Miniconda3-py38_23.5.2-0-Linux-aarch64.sh
```

3. Refresh the shell

```
#Refresh the shell  
source ~/.bashrc  
  
#You should see "(base) username@host" in your terminal.  
#The parenthesis shows you which environment you are currently in.
```

Practice

1. Create a conda environment named training, with python 2.7 as the default python.

```
conda create -n training python=2.7
```

2. Check what other environment is available

```
conda env list
```

3. Activate the environment

```
conda activate training
```

4. Check the python version

```
python --version  
#or  
python  
quit()
```

5. Install R version 3.5

```
conda install -c bioconda rbase=3.5
```

6. Delete the python package

```
conda remove python=2.7
```

7. Verify the uninstallation

```
#try to run python  
python
```

8. Deactivate the environment

```
conda deactivate training
```

9. Delete the environment

```
conda remove -n training -all
```

Prepare conda environment for NGS analysis

1. Create an environment named training with python 3.7 as the main python.

```
conda create -n training python=3.7
```

2. Activate the environment with `conda activate`

3. Install package called micromamba (https://mamba.readthedocs.io/en/latest/user_guide/micromamba.html).

```
conda install -c conda-forge micromamba
```

4. Install the remaining packages using micromamba.

```
micromamba install -c conda-forge -c bioconda \  
blast \  
minimap2 \  
bowtie2 \  
seqkit \  
fastp \  
trimmomatic \  
fastqc \  
bcftools=1.17 \  
samtools=1.17 \  
deepTools \  
nanoplot=1.41 \  
freebayes \  
iVar \  
nanofilt \  
medaka \  
lofreq
```

Information for each package:

blast: <https://www.ncbi.nlm.nih.gov/books/NBK279690/>

minimap2: <https://github.com/lh3/minimap2>

bowtie2: <https://github.com/BenLangmead/bowtie2>

seqkit: <https://github.com/shenwei356/seqkit>
fastp: <https://github.com/OpenGene/fastp>
trimmomatic: <http://www.usadellab.org/cms/?page=trimmomatic>
fastqc: <https://github.com/s-andrews/FastQC>
bcftools: <https://github.com/samtools/bcftools>
samtools: <https://github.com/samtools/samtools>
deepTools: <https://github.com/deeptools/deepTools>
nanoplot: <https://github.com/wdecoster/NanoPlot> [GitHub - wdecoster/NanoPlot: Plotting scripts for long read sequencing data](#)
iVar: <https://andersen-lab.github.io/ivar/html/>
nanofilt: <https://github.com/wdecoster/nanofilt>
medaka: <https://github.com/nanoporetech/medaka>
lofreq: <https://github.com/CSB5/lofreq>